



VARILUX® LENSES AVAILABILITY

for Eyecare Professionals

Effective **JANUARY 2017**



LENSES AVAILABILITY

for Eyecare Professionals

Effective **JANUARY 2017**



| VARILUX COMFORT DRx™ | | LENS ID | | | | | |
|-----------------------------|--|---------|-----------------------------|-----------------|------------------------------------|---------------|----------|
| MATERIAL | PRODUCT | INDEX | COLOR | POWER RANGE | BASE CURVES | ADD POWER | CYLINDER |
| PLASTIC | 1.50 Index Clear | 1.50 | | -10.00 to +6.00 | 2.00, 3.50, 4.25, 6.25, 7.25 | 0.75 to 3.50D | -6.00 |
| | 1.50 Index Essilor Polarized | 1.50 | Gray, Brown, Gray-Green | -10.00 to +6.00 | 2.00, 3.25, 4.25, 6.25, 8.25 | 0.75 to 3.50D | -4.00 |
| | 1.50 Index Xperio UV™ Polarized | 1.50 | Gray, Brown, Gray-Green | -10.00 to +6.00 | 2.00, 3.25, 4.25, 6.25, 8.25 | 0.75 to 3.50D | -4.00 |
| | 1.50 Index Transitions® Signature™ VII | 1.50 | Graphite Green, Gray, Brown | -10.00 to +6.00 | 2.00, 3.50, 4.25, 6.25, 7.25 | 0.75 to 3.50D | -6.00 |
| | 1.50 Index Transitions® XTRActive™ | 1.50 | Gray | -10.00 to +6.00 | 2.25, 4.25, 5.25, 8.25 | 0.75 to 3.50D | -6.00 |
| | 1.50 Index Transitions® Vantage™ | 1.50 | Gray | -10.00 to +6.00 | 2.25, 4.25, 5.25, 8.25 | 0.75 to 3.50D | -6.00 |
| TREXA™ (TRIVEX) | Clear | 1.53 | | -10.00 to +6.00 | 2.00, 3.50, 4.25, 6.25, 7.25 | 0.75 to 3.50D | -6.00 |
| | Transitions Signature VII | 1.53 | Gray, Brown | -10.00 to +6.00 | 2.00, 3.50, 4.25, 6.25, 7.25 | 0.75 to 3.50D | -6.00 |
| | Transitions XTRActive | 1.53 | Graphite Green, Gray, Brown | -10.00 to +6.00 | 2.00, 3.50, 4.25, 6.25, 7.25 | 0.75 to 3.50D | -6.00 |
| AIRWEAR™ (POLYCARBONATE) | Clear | 1.59 | | -10.00 to +8.00 | 1.50, 3.25, 4.25, 5.50, 7.00, 8.50 | 0.75 to 3.50D | -6.00 |
| | Essilor Polarized | 1.59 | Gray, Brown, Gray-Green | -10.00 to +8.00 | 1.50, 3.25, 4.25, 5.50, 7.00, 8.50 | 0.75 to 3.50D | -6.00 |
| | Xperio UV Polarized | 1.59 | Gray, Brown, Gray-Green | -10.00 to +8.00 | 1.50, 3.25, 4.25, 5.50, 7.00, 8.50 | 0.75 to 3.50D | -6.00 |
| | Transitions Signature VII | 1.59 | Graphite Green, Gray, Brown | -10.00 to +8.00 | 1.50, 3.25, 4.25, 5.50, 7.00, 8.50 | 0.75 to 3.50D | -6.00 |
| | Transitions XTRActive | 1.59 | Gray | -10.00 to +8.00 | 1.50, 3.25, 4.25, 5.50, 7.00, 8.50 | 0.75 to 3.50D | -6.00 |
| | Transitions Vantage | 1.59 | Gray | -10.00 to +8.00 | 1.50, 3.25, 4.25, 5.50, 7.00, 8.50 | 0.75 to 3.50D | -6.00 |
| THIN&LITE™ | 1.60 Clear | 1.60 | | -12.00 to +8.00 | 1.75, 2.75, 3.75, 6.00, 6.75, 8.50 | 0.75 to 3.50D | -6.00 |
| | 1.60 Transitions Signature VII | 1.60 | Gray, Brown | -12.00 to +8.00 | 1.75, 2.75, 3.75, 6.00, 6.75, 8.50 | 0.75 to 3.50D | -6.00 |
| | 1.67 Clear | 1.67 | | -12.00 to +9.00 | 1.75, 2.75, 3.75, 6.00, 6.75, 8.00 | 0.75 to 3.50D | -6.00 |
| | 1.67 Essilor Polarized | 1.67 | Gray, Brown, Gray-Green | -12.00 to +9.00 | 1.75, 2.75, 3.75, 5.25, 6.75, 8.00 | 0.75 to 3.50D | -6.00 |
| | 1.67 Xperio UV Polarized | 1.67 | Gray, Brown, Gray-Green | -12.00 to +9.00 | 1.75, 2.75, 3.75, 5.25, 6.75, 8.00 | 0.75 to 3.50D | -4.00 |
| | 1.67 Transitions Signature VII | 1.67 | Graphite Green, Gray, Brown | -12.00 to +9.00 | 1.75, 2.75, 3.75, 6.00, 6.75, 8.00 | 0.75 to 3.50D | -4.00 |
| | 1.67 Transitions XTRActive | 1.67 | Gray | -12.00 to +9.00 | 1.75, 2.75, 3.75, 6.00, 6.75, 8.00 | 0.75 to 3.50D | -6.00 |

MFH=17mm

NOTES: All Varilux digital lenses require a hard coat. Essilor AR stacks are required when No-Glare coating is applied (for example, Crizal® No-Glare lenses). Smart Blue Filter™ feature will be included in clear TREXA™ (TRIVEX) and clear Thin&Lite™ 1.74 in February 2017.

| VARILUX COMFORT SHORT DRx™ | | LENS ID | | | | | |
|--|---------|---------|-------|-------------|-------------|-----------|----------|
| MATERIAL | PRODUCT | INDEX | COLOR | POWER RANGE | BASE CURVES | ADD POWER | CYLINDER |
| Rx availability for Varilux Comfort Short DRx matches Varilux Comfort DRx. | | | | | | | |

MFH=14mm

VARILUX.COM



| VARILUX S DESIGN™ | | LENS ID | | VARILUX S FIT™ | | LENS ID | | VARILUX S 4D™ | | LENS ID | |
|-----------------------------|--|---------|-----------------------------|------------------|---------------|--------------|--|---------------|--|---------|--|
| MATERIAL | PRODUCT | INDEX | COLOR | POWER RANGE | ADD POWER* | MAX CYLINDER | | | | | |
| PLASTIC | 1.50 Index Clear | 1.50 | | -10.00 to +6.00 | 0.75 to 4.00D | -6.00 | | | | | |
| | 1.50 Index Essilor Polarized | 1.50 | Gray, Brown | -10.00 to +6.00 | 0.75 to 4.00D | -6.00 | | | | | |
| | 1.50 Index Xperio UV™ Polarized | 1.50 | Gray, Brown | -10.00 to +6.00 | 0.75 to 4.00D | -6.00 | | | | | |
| | 1.50 Index Transitions® Signature™ VII | 1.50 | Gray, Brown | -10.00 to +6.00 | 0.75 to 4.00D | -6.00 | | | | | |
| | 1.50 Index Transitions® XTRActive™ | 1.50 | Gray, Brown | -10.00 to +6.00 | 0.75 to 4.00D | -6.00 | | | | | |
| TREXA™ (TRIVEX) | Clear | 1.53 | | -10.00 to +6.00 | 0.75 to 4.00D | -6.00 | | | | | |
| | Transitions Signature VII | 1.53 | Gray | -10.00 to +6.00 | 0.75 to 4.00D | -6.00 | | | | | |
| AIRWEAR™ (POLYCARBONATE) | Clear | 1.59 | | -10.00 to +6.00 | 0.75 to 4.00D | -6.00 | | | | | |
| | Essilor Polarized | 1.59 | Gray, Brown | -10.00 to +6.00 | 0.75 to 4.00D | -6.00 | | | | | |
| | Xperio UV Polarized | 1.59 | Gray, Brown | -10.00 to +6.00 | 0.75 to 4.00D | -6.00 | | | | | |
| | Transitions Signature VII | 1.59 | Graphite Green, Gray, Brown | -10.00 to +6.00 | 0.75 to 4.00D | -6.00 | | | | | |
| | Transitions XTRActive | 1.59 | Graphite Green, Gray, Brown | -10.00 to +6.00 | 0.75 to 4.00D | -6.00 | | | | | |
| THIN&LITE™ | 1.67 Clear | 1.67 | | -12.00 to +9.00 | 0.75 to 4.00D | -6.00 | | | | | |
| | 1.67 Transitions Signature VII | 1.67 | Graphite Green, Gray, Brown | -12.00 to +9.00 | 0.75 to 4.00D | -6.00 | | | | | |
| | 1.67 Transitions XTRActive | 1.67 | Graphite Green, Gray, Brown | -12.00 to +9.00 | 0.75 to 4.00D | -6.00 | | | | | |
| | 1.74 Clear** | 1.74 | | -14.00 to +13.00 | 0.75 to 4.00D | -6.00 | | | | | |

MFH=14mm

NOTES: *ADD Power in steps of 0.25D. At fitting height less than 17mm, ADD power is up to +3.50. Thin&Lite 1.74 lenses come systematic with choice of Crizal® Previncia®, Crizal Sapphire UV™, Crizal Avancé UV™ or Crizal Alizé UV™. All Varilux digital lenses require a hard coat. Essilor AR stacks are required when No-Glare coating is applied (for example, Crizal® No-Glare lenses).

** Varilux S Series available in Thin&Lite 1.74 Clear with Smart Blue Filter feature in February 2017.

VARILUX® LENSES AVAILABILITY

for Eyecare Professionals

Effective **JANUARY 2017**

VARILUX® LENSES AVAILABILITY

for Eyecare Professionals

Effective **JANUARY 2017**

| VARILUX® PHYSIO® W3+ | | TEMPORAL | NASAL | VARILUX® PHYSIO® W3+ FIT | | TEMPORAL | NASAL | VARILUX® PHYSIO® W3+ EYECODE™ | | TEMPORAL | NASAL |
|--------------------------------|--|----------|-----------------------------|---|---|---------------------|----------|-------------------------------|--|----------|-------|
| MATERIAL | PRODUCT | INDEX | COLOR | POWER RANGE | BASE CURVES | ADD POWER* | CYLINDER | | | | |
| PLASTIC | 1.50 Index Clear | 1.50 | | -10.00 to +6.00 | 2.00, 3.25, 4.50, 5.75, 7.25 | 0.75 to 3.50, 4.00D | -6.00 | | | | |
| | 1.50 Index Essilor Polarized | 1.50 | Gray, Brown | -10.00 to +6.00 | 2.00, 3.25, 4.50, 5.75, 7.25 | 0.75 to 3.50D | -6.00 | | | | |
| | 1.50 Index Xperio UV™ Polarized | 1.50 | Gray, Brown | -10.00 to +6.00 | 2.00, 3.25, 4.50, 5.75, 7.25 | 0.75 to 3.50D | -6.00 | | | | |
| | 1.50 Index Transitions® Signature™ VII | 1.50 | Gray, Brown | -10.00 to +6.00 | 2.00, 3.25, 4.50, 5.75, 7.25 | 0.75 to 3.50, 4.00D | -6.00 | | | | |
| | 1.50 Index Transitions® XTRActive® | 1.50 | Gray | -10.00 to +6.00 | 2.00, 3.25, 4.50, 5.75, 7.25 | 0.75 to 3.50, 4.00D | -6.00 | | | | |
| TREXA™ (TRIVEX™) | Clear | 1.53 | | -10.00 to +6.00 | 2.00, 3.25, 4.50, 5.75, 7.25 | 0.75 to 3.50, 4.00D | -6.00 | | | | |
| | Transitions Signature VII | 1.53 | Gray, Brown | -10.00 to +6.00 | 2.00, 3.25, 4.50, 5.75, 7.25 | 0.75 to 3.50, 4.00D | -6.00 | | | | |
| AIRWEAR (POLYCARBONATE) | Clear | 1.59 | | -10.00 to +6.00 | 1.75, 3.00, 4.00, 5.50, 7.25 | 0.75 to 3.50, 4.00D | -6.00 | | | | |
| | Essilor Polarized | 1.59 | Gray, Brown | -10.00 to +6.00 | 1.75, 3.00, 4.00, 5.50, 7.25 | 0.75 to 3.50D | -6.00 | | | | |
| | Xperio UV Polarized | 1.59 | Gray, Brown | -10.00 to +6.00 | 1.75, 3.00, 4.00, 5.50, 7.25 | 0.75 to 3.50D | -6.00 | | | | |
| | Transitions Signature VII | 1.59 | Graphite Green, Gray, Brown | -10.00 to +6.00 | 1.75, 3.00, 4.00, 5.50, 7.25 | 0.75 to 3.50, 4.00D | -6.00 | | | | |
| | Transitions XTRActive | 1.59 | Gray, Brown | -10.00 to +6.00 | 1.75, 3.00, 4.00, 5.50, 7.25 | 0.75 to 3.50, 4.00D | -6.00 | | | | |
| THIN&LITE | 1.60 Clear | 1.60 | | -12.00 to +8.00 | 1.75, 3.00, 4.00, 5.50, 7.00, 8.00 | 0.75 to 3.50, 4.00D | -6.00 | | | | |
| | 1.60 Transitions Signature VII | 1.60 | Brown | -12.00 to +8.00 | 1.75, 3.00, 4.00, 5.50, 7.00, 8.00 | 0.75 to 3.50, 4.00D | -6.00 | | | | |
| | 1.67 Clear | 1.67 | | -14.00 to +9.00 | 1.75, 2.75, 3.75, 5.25, 6.50, 8.00 | 0.75 to 3.50, 4.00D | -6.00 | | | | |
| | 1.67 Transitions Signature VII | 1.67 | Gray, Brown | -14.00 to +9.00 | 1.75, 2.75, 3.75, 5.25, 6.50, 8.00 | 0.75 to 3.50, 4.00D | -6.00 | | | | |
| | 1.67 Transitions XTRActive | 1.67 | Gray | -12.00 to +9.00 | 1.75, 2.75, 3.75, 5.25, 6.50, 8.00 | 0.75 to 3.50, 4.00D | -6.00 | | | | |
| | 1.74 Clear | 1.74 | | -18.00 to +13.00 | 1.75, 2.75, 3.75, 5.25, 6.50, 8.00, 10.00 | 0.75 to 3.50, 4.00D | -6.00 | | | | |
| 1.74 Transitions Signature VII | 1.74 | Gray | -18.00 to +13.00 | 1.75, 2.75, 3.75, 5.25, 6.50, 8.00, 10.00 | 0.75 to 3.50, 4.00D | -6.00 | | | | | |

MFH=14mm

| VARILUX COMFORT® W2+ | | TEMPORAL | NASAL | VARILUX COMFORT® W2+ FIT | | TEMPORAL | NASAL | VARILUX COMFORT® W2+ EYECODE™ | | TEMPORAL | NASAL |
|----------------------------|--------------------------------|----------|-------------|--------------------------|------------------------------------|---------------|----------|-------------------------------|--|----------|-------|
| MATERIAL | PRODUCT | INDEX | COLOR | POWER RANGE | BASE CURVES | ADD POWER* | CYLINDER | | | | |
| PLASTIC | 1.50 Index Clear | 1.50 | | -10.00 to +6.00 | 2.00, 3.25, 4.50, 5.75, 7.25 | 0.75 to 3.50D | -6.00 | | | | |
| | 1.50 Transitions Signature VII | 1.50 | Gray, Brown | -10.00 to +6.00 | 2.00, 3.25, 4.50, 5.75, 7.25 | 0.75 to 3.50D | -6.00 | | | | |
| | 1.50 Transitions XTRActive | 1.50 | Gray | -10.00 to +6.00 | 2.00, 3.25, 4.50, 5.75, 7.25 | 0.75 to 3.50D | -6.00 | | | | |
| TREXA™ (TRIVEX™) | Clear | 1.53 | | -10.00 to +6.00 | 2.00, 3.25, 4.50, 5.75, 7.25 | 0.75 to 3.50D | -6.00 | | | | |
| | Transitions Signature VII | 1.53 | Gray, Brown | -10.00 to +6.00 | 2.00, 3.25, 4.50, 5.75, 7.25 | 0.75 to 3.50D | -6.00 | | | | |
| AIRWEAR (POLYCARBONATE) | Clear | 1.59 | | -10.00 to +6.00 | 1.75, 3.00, 4.00, 5.50, 7.00 | 0.75 to 3.50D | -6.00 | | | | |
| | Essilor Polarized | 1.59 | Gray, Brown | -10.00 to +6.00 | 1.75, 3.00, 4.00, 5.50, 7.00 | 0.75 to 3.50D | -6.00 | | | | |
| | Xperio UV Polarized | 1.59 | Gray, Brown | -10.00 to +6.00 | 1.75, 3.00, 4.00, 5.50, 7.00 | 0.75 to 3.50D | -6.00 | | | | |
| | Transitions Signature VII | 1.59 | Gray, Brown | -10.00 to +6.00 | 1.75, 3.00, 4.00, 5.50, 7.00 | 0.75 to 3.50D | -6.00 | | | | |
| | Transitions XTRActive | 1.59 | Gray, Brown | -10.00 to +6.00 | 1.75, 3.00, 4.00, 5.50, 7.00 | 0.75 to 3.50D | -6.00 | | | | |
| THIN&LITE | 1.60 Clear | 1.60 | | -12.00 to +8.00 | 1.75, 3.00, 4.00, 5.50, 7.00, 8.00 | 0.75 to 3.50D | -6.00 | | | | |
| | 1.60 Transitions Signature VII | 1.60 | Brown | -12.00 to +8.00 | 1.75, 3.00, 4.00, 5.50, 7.00, 8.00 | 0.75 to 3.50D | -6.00 | | | | |
| | 1.67 Clear | 1.67 | | -14.00 to +9.00 | 1.75, 2.75, 3.75, 5.25, 6.50, 8.00 | 0.75 to 3.50D | -6.00 | | | | |
| | 1.67 Transitions Signature VII | 1.67 | Gray, Brown | -14.00 to +9.00 | 1.75, 2.75, 3.75, 5.25, 6.50, 8.00 | 0.75 to 3.50D | -6.00 | | | | |
| | 1.67 Transitions XTRActive | 1.67 | Gray | -12.00 to +9.00 | 1.75, 2.75, 3.75, 5.25, 6.50, 8.00 | 0.75 to 3.00D | -6.00 | | | | |

MFH=14mm

NOTES: *ADD Power in steps of 0.25D. At fitting height less than 17mm, ADD power is up to +3.00.

Thin&Lite 1.74 lenses come systematic with choice of Crizal® Previncia®, Crizal Sapphire UV™, Crizal Avancé UV™ or Crizal Alizé UV™

All Varilux digital lenses require a hard coat. Essilor AR stacks are required when No-Glare coating is applied (for example, Crizal® No-Glare lenses).

Varilux Comfort W2+ and Varilux Physio W3+ lenses in clear TREXA™ (TRIVEX™), clear Thin&Lite 1.60 will have Smart Blue Filter™ feature in February 2017.

**TREXA™ (TRIVEX™) material will be available with Physio W3+ beginning in February 2017.

| VARILUX PHYSIO DRX™ | | LENS ID | VARILUX PHYSIO DRX™ | | LENS ID | | |
|----------------------------|--|---------|-----------------------------|------------------|---|---------------------|----------|
| MATERIAL | PRODUCT | INDEX | COLOR | POWER RANGE | BASE CURVES | ADD POWER | CYLINDER |
| PLASTIC | 1.50 Index Clear | 1.50 | | -10.00 to +6.00 | 2.00, 3.50, 4.25, 6.25, 7.25 | 0.75 to 3.50, 4.00D | -6.00 |
| | 1.50 Index Essilor Polarized | 1.50 | Gray, Brown, Gray-Green | -10.00 to +6.00 | 2.00, 3.25, 4.25, 6.25, 8.25 | 0.75 to 3.50D | -4.00 |
| | 1.50 Index Xperio UV™ Polarized | 1.50 | Gray, Brown, Gray-Green | -10.00 to +6.00 | 2.00, 3.25, 4.25, 6.25, 8.25 | 0.75 to 3.50D | -4.00 |
| | 1.50 Index Transitions® Signature™ VII | 1.50 | Graphite Green, Gray, Brown | -10.00 to +6.00 | 2.00, 3.50, 4.25, 6.25, 7.25 | 0.75 to 3.50, 4.00D | -6.00 |
| | 1.50 Index Transitions® XTRActive® | 1.50 | Gray | -10.00 to +6.00 | 2.25, 4.25, 5.25, 8.25 | 0.75 to 3.50, 4.00D | -6.00 |
| | 1.50 Index Transitions® Vantage™ | 1.50 | Gray | -10.00 to +6.00 | 2.25, 4.25, 5.25, 8.25 | 0.75 to 3.50, 4.00D | -6.00 |
| TREXA™ (TRIVEX™) | Clear | 1.53 | | -10.00 to +6.00 | 2.00, 3.50, 4.25, 6.25, 7.25 | 0.75 to 3.50, 4.00D | -6.00 |
| | Transitions Signature VII | 1.53 | Gray, Brown | -10.00 to +6.00 | 2.00, 3.50, 4.25, 6.25, 7.25 | 0.75 to 3.50, 4.00D | -6.00 |
| AIRWEAR (POLYCARBONATE) | Clear | 1.59 | | -10.00 to +8.00 | 1.50, 3.25, 4.25, 5.50, 7.00, 8.50 | 0.75 to 3.50, 4.00D | -6.00 |
| | Essilor Polarized | 1.59 | Gray, Brown, Gray-Green | -10.00 to +8.00 | 1.50, 3.25, 4.25, 5.50, 7.00, 8.50 | 0.75 to 3.50, 4.00D | -6.00 |
| | Xperio UV Polarized | 1.59 | Gray, Brown, Gray-Green | -10.00 to +8.00 | 1.50, 3.25, 4.25, 5.50, 7.00, 8.50 | 0.75 to 3.50, 4.00D | -6.00 |
| | Transitions Signature VII | 1.59 | Graphite Green, Gray, Brown | -10.00 to +8.00 | 1.50, 3.25, 4.25, 5.50, 7.00, 8.50 | 0.75 to 3.50, 4.00D | -6.00 |
| | Transitions XTRActive | 1.59 | Gray | -10.00 to +8.00 | 1.50, 3.25, 4.25, 5.50, 7.00, 8.50 | 0.75 to 3.50, 4.00D | -6.00 |
| THIN&LITE | 1.60 Clear | 1.60 | | -12.00 to +8.00 | 1.75, 2.75, 3.75, 6.00, 6.75, 8.50 | 0.75 to 3.50, 4.00D | -6.00 |
| | 1.60 Transitions Signature VII | 1.60 | Gray, Brown | -12.00 to +8.00 | 1.75, 2.75, 3.75, 6.00, 6.75, 8.50 | 0.75 to 3.50, 4.00D | -6.00 |
| | 1.67 Clear | 1.67 | | -12.00 to +9.00 | 1.75, 2.75, 3.75, 6.00, 6.75, 8.00 | 0.75 to 3.50, 4.00D | -6.00 |
| | 1.67 Essilor Polarized | 1.67 | Gray, Brown, Gray-Green | -12.00 to +9.00 | 1.75, 2.75, 3.75, 5.25, 6.75, 8.00 | 0.75 to 3.50, 4.00D | -6.00 |
| | 1.67 Xperio UV Polarized | 1.67 | Gray, Brown, Gray-Green | -12.00 to +9.00 | 1.75, 2.75, 3.75, 5.25, 6.75, 8.00 | 0.75 to 3.50, 4.00D | -6.00 |
| | 1.67 Transitions Signature VII | 1.67 | Graphite Green, Gray, Brown | -12.00 to +9.00 | 1.75, 2.75, 3.75, 6.00, 6.75, 8.00 | 0.75 to 3.50, 4.00D | -6.00 |
| | 1.67 Transitions XTRActive | 1.67 | Gray | -12.00 to +9.00 | 1.75, 2.75, 3.75, 6.00, 6.75, 8.00 | 0.75 to 3.50, 4.00D | -6.00 |
| | 1.74 Clear | 1.74 | | -20.00 to +13.00 | 0.75, 1.75, 2.75, 3.75, 5.25, 6.50, 8.00, 10.00 | 0.75 to 3.50, 4.00D | -6.00 |
| | 1.74 Transitions Signature VII | 1.74 | Gray | -20.00 to +13.00 | 0.75, 1.75, 2.75, 3.75, 5.25, 6.50, 8.00, 10.00 | 0.75 to 3.50, 4.00D | -6.00 |

MFH=17mm

NOTES: *Thin&Lite 1.74 lenses come systematic with choice of Crizal® Previncia®, Crizal Sapphire UV™, Crizal Avancé UV™ or Crizal Alizé UV™
Varilux Physio DRx lenses in clear TREXA™ (TRIVEX™), and clear Thin&Lite 1.74 will have Smart Blue Filter™ feature later in 2016.
Smart Blue Filter™ feature will be included in clear TREXA™ (TRIVEX™) and clear Thin&Lite 1.74 in February 2017.

VARILUX PHYSIO SHORT DRX™

Rx availability for Varilux Physio Short DRx matches Varilux Physio DRx, except ADD power range is 0.75 to 3.00D across all materials.

MFH=14mm